

**ELITE VETERINARY PRODUCT'S
GAME READY CANINE SPRING DEAL!**

Complete Canine system

\$3099.00

List price \$3,750.00

Helping dogs heal faster with cold and compression

Ask us about our Post Surgical rental program, now your clients can have the same state of the art therapy, at home!

*Elite Veterinary Products
Your source for Game Ready*

Elite Veterinary Products

7828 Percussion Dr

Apex, NC 27539

919-800-9044

jdonovin@elitevetproducts.com

www.elitevetproducts.com

**Elite Vet Products is an authorized distributor
of the Game Ready System for both Canine and Equine systems.**

G A M E  R E A D Y ®

What is Game ready and what does it do?

Game Ready is proven therapy - Based upon a parallel system for humans, Game Ready combines: Dry cold therapy and active compression in one revolutionary application.

Cold therapy and minimize pain and muscle spasm, reduce swelling and decrease damage to healthy tissue following an injury or surgery.

Compression can reduce swelling, prevent tissue damage and promote faster healing.

In fact, Game Ready is so effective it is already the injury and treatment system by choice of pro athletes, collegiate sports teams and the leading orthopedic clinics. It is also used by the top canine and equine veterinarians around the world.

Discover a non-invasive way to help animals heal faster from injury or surgery.

How does it work?

Game Ready delivers a unique combination of adjustable cold therapy and active compression via a microprocessor-regulated control unit and formfitting dual action wraps.



Follow the steps below to use the Game Ready System

- Fill unit with 3/4 ice and cold water
- Plug unit into AC outlet
- Place appropriate wrap around injured area. The hose connector **MUST** be facing up. Ensure that the wrap is secured snug around the area.
- Connect the hose from the wrap to the control unit. The hose only connects one way, you will hear a click when connected properly.
- Press the POWER button
- Set PRESSURE to medium (do not use High pressure, if dog is sensitive to the pressure, change to low). Pressure button is first left button on top
- Set RUN TIME to 25 minutes - Run Time is the second button on the bottom
- Press the PLAY button
- Use system at least 3x times a day or as recommended by your veterinarian

Replace ice and cold water after each use.